



Protein PR73 Antibody

Product Code	CSB-PA361020XA01MIV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03320
Immunogen	Recombinant Mouse mammary tumor virus (strain GR) (MMTV) Protein PR73 protein
Raised In	Rabbit
Species Reactivity	Mouse mammary tumor virus (strain GR) (MMTV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mouse mammary tumor virus (strain GR) (MMTV)
Lead Time	Made-to-order (14-16 weeks)