



sag Antibody

Product Code	CSB-PA365906XA01MIU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03319
Immunogen	Recombinant Mouse mammary tumor virus (strain C3H) (MMTV) sag protein
Raised In	Rabbit
Species Reactivity	Mouse mammary tumor virus (strain C3H) (MMTV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mouse mammary tumor virus (strain C3H) (MMTV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	sag
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.