



RAS Antibody

Product Code	CSB-PA365522XA01RBR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01114
Immunogen	Recombinant Rasheed rat sarcoma virus RAS protein
Raised In	Rabbit
Species Reactivity	Rasheed rat sarcoma virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rasheed rat sarcoma virus
Lead Time	Made-to-order (14-16 weeks)
Target Names	RAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.