



rev Antibody

Product Code	CSB-PA335653XA01CDE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31628
Immunogen	Recombinant Caprine arthritis encephalitis virus (CAEV) rev protein
Raised In	Rabbit
Species Reactivity	Caprine arthritis encephalitis virus (CAEV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Caprine arthritis encephalitis virus (CAEV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rev
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.