



rev Antibody

Product Code	CSB-PA332964XA01BKB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24097
Immunogen	Recombinant Bovine immunodeficiency virus (strain R29) (BIV) (Bovine immunodeficiency-like virus) rev protein
Raised In	Rabbit
Species Reactivity	Bovine immunodeficiency virus (strain R29) (BIV) (Bovine immunodeficiency-like virus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bovine immunodeficiency virus (strain R29) (BIV) (Bovine immunodeficiency-like virus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rev
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.