



rev Antibody

Product Code	CSB-PA324983XA01FAK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20885
Immunogen	Recombinant Feline immunodeficiency virus (isolate Petaluma) (FIV) rev protein
Raised In	Rabbit
Species Reactivity	Feline immunodeficiency virus (isolate Petaluma) (FIV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Feline immunodeficiency virus (isolate Petaluma) (FIV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rev
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.