



V-SRC Antibody

Product Code	CSB-PA365441XA01RHO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00526
Immunogen	Recombinant Rous sarcoma virus (strain Prague C) (RSV-PrC) V-SRC protein
Raised In	Rabbit
Species Reactivity	Rous sarcoma virus (strain Prague C) (RSV-PrC)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rous sarcoma virus (strain Prague C) (RSV-PrC)
Lead Time	Made-to-order (14-16 weeks)
Target Names	V-SRC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.