



vif Antibody

Product Code	CSB-PA767623XA01FYS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q82852
Immunogen	Recombinant Jembrana disease virus (JDV) vif protein
Raised In	Rabbit
Species Reactivity	Jembrana disease virus (JDV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Jembrana disease virus (JDV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	vif
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.