



vpu Antibody

Product Code	CSB-PA327066XA01SFU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31624
Immunogen	Recombinant Sheep pulmonary adenomatosis virus (Jaagsiekte sheep retrovirus) (JSRV) vpu protein
Raised In	Rabbit
Species Reactivity	Sheep pulmonary adenomatosis virus (Jaagsiekte sheep retrovirus) (JSRV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sheep pulmonary adenomatosis virus (Jaagsiekte sheep retrovirus) (JSRV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	vpu
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.