



vif Antibody

Product Code	CSB-PA301699XA01MAW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69715
Immunogen	Recombinant Maedi visna virus (strain 1514) (MVV) (Visna lentivirus) vif protein
Raised In	Rabbit
Species Reactivity	Maedi visna virus (strain 1514) (MVV) (Visna lentivirus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Maedi visna virus (strain 1514) (MVV) (Visna lentivirus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	vif
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.