



Uncharacterized protein ORF3 Antibody

Product Code	CSB-PA321658XA01FAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19031
Immunogen	Recombinant Feline immunodeficiency virus (strain San Diego) (FIV) Uncharacterized protein ORF3 protein
Raised In	Rabbit
Species Reactivity	Feline immunodeficiency virus (strain San Diego) (FIV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Feline immunodeficiency virus (strain San Diego) (FIV)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.