



# HSP104 Antibody

|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA329079XA01SVG                                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                    |
| Uniprot No.         | P31539                                                                                           |
| Immunogen           | Recombinant Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) HSP104 protein |
| Raised In           | Rabbit                                                                                           |
| Species Reactivity  | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                            |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                                     |
| Form                | Liquid                                                                                           |
| Conjugate           | Non-conjugated                                                                                   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                 |
| Purification Method | Antigen Affinity Purified                                                                        |
| Isotype             | IgG                                                                                              |
| Clonality           | Polyclonal                                                                                       |
| Product Type        | Polyclonal Antibody                                                                              |
| Immunogen Species   | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)                            |
| Lead Time           | Made-to-order (14-16 weeks)                                                                      |
| Target Names        | HSP104                                                                                           |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                      |