



abf1 Antibody

Product Code	CSB-PA835663XA01HYE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92455
Immunogen	Recombinant Hypocrea jecorina (Trichoderma reesei) abf1 protein
Raised In	Rabbit
Species Reactivity	Hypocrea jecorina (Trichoderma reesei)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Hypocrea jecorina (Trichoderma reesei)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	abf1