



ADH8 Antibody

Product Code	CSB-PA530020XA01ETD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O57380
Immunogen	Recombinant Pelophylax perezi (Perez's frog) (Rana perezi) ADH8 protein
Raised In	Rabbit
Species Reactivity	Pelophylax perezi (Perez's frog) (Rana perezi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pelophylax perezi (Perez's frog) (Rana perezi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ADH8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.