



LKA1 Antibody

Product Code	CSB-PA312441XA01LGO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01117
Immunogen	Recombinant Lipomyces kononenkoae (Yeast) LKA1 protein
Raised In	Rabbit
Species Reactivity	Lipomyces kononenkoae (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lipomyces kononenkoae (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	LKA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.