



GAMP Antibody

Product Code	CSB-PA291826XA01FZO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03045
Immunogen	Recombinant Amorphytheca resiniae (Creosote fungus) (Hormoconis resiniae) GAMP protein
Raised In	Rabbit
Species Reactivity	Amorphytheca resiniae (Creosote fungus) (Hormoconis resiniae)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Amorphytheca resiniae (Creosote fungus) (Hormoconis resiniae)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GAMP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.