



GAA Antibody

Product Code	CSB-PA336965XA01BTJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42042
Immunogen	Recombinant Blastobotrys adenivorans (Yeast) (Arxula adenivorans) GAA protein
Raised In	Rabbit
Species Reactivity	Blastobotrys adenivorans (Yeast) (Arxula adenivorans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Blastobotrys adenivorans (Yeast) (Arxula adenivorans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GAA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.