



amy Antibody

Product Code	CSB-PA179916XA01APS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30292
Immunogen	Recombinant Aspergillus shirousami amy protein
Raised In	Rabbit
Species Reactivity	Aspergillus shirousami
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aspergillus shirousami
Lead Time	Made-to-order (14-16 weeks)
Gene Names	amy