



B1 Antibody

Product Code	CSB-PA303087XA01FCA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68829
Immunogen	Recombinant Flock house virus (FHV) B1 protein
Raised In	Rabbit
Species Reactivity	Flock house virus (FHV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Flock house virus (FHV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.