



BBP1 Antibody

Product Code	CSB-PA520505XA01SVJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7QA20
Immunogen	Recombinant Saccharomyces cerevisiae (strain FostersB) (Baker's yeast) BBP1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain FostersB) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain FostersB) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	BBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.