



CDC25 Antibody

Product Code	CSB-PA589917XA01LMA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02342
Immunogen	Recombinant Lachancea kluyveri (Yeast) (Saccharomyces kluyveri) CDC25 protein
Raised In	Rabbit
Species Reactivity	Lachancea kluyveri (Yeast) (Saccharomyces kluyveri)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lachancea kluyveri (Yeast) (Saccharomyces kluyveri)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CDC25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.