



CYP52A4 Antibody

Product Code	CSB-PA153521XA01CZJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16141
Immunogen	Recombinant Candida maltosa (Yeast) CYP52A4 protein
Raised In	Rabbit
Species Reactivity	Candida maltosa (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida maltosa (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CYP52A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.