



CRP1 Antibody

Product Code	CSB-PA520684XA01SVO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7KPJ0
Immunogen	Recombinant Saccharomyces cerevisiae (strain Lalvin QA23) (Baker's yeast) CRP1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain Lalvin QA23) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain Lalvin QA23) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.