



Cytochrome b5 Antibody

Product Code	CSB-PA862097XA01RDF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HFV1
Immunogen	Recombinant Rhizopus stolonifer (Rhizopus nigricans) Cytochrome b5 protein
Raised In	Rabbit
Species Reactivity	Rhizopus stolonifer (Rhizopus nigricans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhizopus stolonifer (Rhizopus nigricans)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.