



Cyp2d4 Antibody

Product Code	CSB-PA714505XA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q64680
Immunogen	Recombinant Rattus norvegicus (Rat) Cyp2d4 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus (Rat)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Cyp2d4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.