



# Fad5 Antibody

Product Code	CSB-PA443402XA01FZU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8S3C1
Immunogen	Recombinant Thraustochytrium sp. Fad5 protein
Raised In	Rabbit
Species Reactivity	Thraustochytrium sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thraustochytrium sp.
Lead Time	Made-to-order (14-16 weeks)
Gene Names	Fad5