



Fad4 Antibody

Product Code	CSB-PA824892XA01FZU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8S3C0
Immunogen	Recombinant Thraustochytrium sp. Fad4 protein
Raised In	Rabbit
Species Reactivity	Thraustochytrium sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thraustochytrium sp.
Lead Time	Made-to-order (14-16 weeks)
Target Names	Fad4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.