



POL2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA018309XA01SVG |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P21951 |
| Immunogen | Synthetic Peptide (1320-1344aa) |
| Species Reactivity | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Affinity-chromatography |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) |
| Target Names | POL2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |