



SUP2 Antibody

Product Code	CSB-PA340559XA01OAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23637
Immunogen	Recombinant Ogataea pini (Yeast) (Pichia pinus) SUP2 protein
Raised In	Rabbit
Species Reactivity	Ogataea pini (Yeast) (Pichia pinus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ogataea pini (Yeast) (Pichia pinus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	SUP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.