



Formate dehydrogenase Antibody

Product Code	CSB-PA333773XA01PVB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33160
Immunogen	Recombinant Pseudomonas sp. (strain 101) (Achromobacter parvulus T1) Formate dehydrogenase protein
Raised In	Rabbit
Species Reactivity	Pseudomonas sp. (strain 101) (Achromobacter parvulus T1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas sp. (strain 101) (Achromobacter parvulus T1)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.