



Delta(12) fatty acid desaturase Antibody

Product Code	CSB-PA896802XA01FTD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y8H5
Immunogen	Recombinant Mortierella alpina (Oleaginous fungus) (Mortierella renispora) Delta(12) fatty acid desaturase protein
Raised In	Rabbit
Species Reactivity	Mortierella alpina (Oleaginous fungus) (Mortierella renispora)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mortierella alpina (Oleaginous fungus) (Mortierella renispora)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.