



MET7 Antibody

Product Code	CSB-PA008851XA01SVR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7QKX4
Immunogen	Recombinant Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast) MET7 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	MET7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.