



GPA1 Antibody

Product Code	CSB-PA517319XA01EWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O04278
Immunogen	Recombinant Pisum sativum (Garden pea) GPA1 protein
Raised In	Rabbit
Species Reactivity	Pisum sativum (Garden pea)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pisum sativum (Garden pea)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GPA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.