



grisea Antibody

Product Code	CSB-PA835639XA01EXI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92258
Immunogen	Recombinant Podospora anserina (Pleurance anserina) grisea protein
Raised In	Rabbit
Species Reactivity	Podospora anserina (Pleurance anserina)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Podospora anserina (Pleurance anserina)
Lead Time	Made-to-order (14-16 weeks)
Target Names	grisea
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.