



cekA Antibody

Product Code	CSB-PA613252XA01APO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12679
Immunogen	Recombinant Aspergillus kawachii (strain NBRC 4308) (White koji mold) (Aspergillus awamori var. kawachi) cekA protein
Raised In	Rabbit
Species Reactivity	Aspergillus kawachii (strain NBRC 4308) (White koji mold) (Aspergillus awamori var. kawachi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aspergillus kawachii (strain NBRC 4308) (White koji mold) (Aspergillus awamori var. kawachi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	cekA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.