



ORM1 Antibody

Product Code	CSB-PA856867XA01FZM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96495
Immunogen	Recombinant Saccharomyces pastorianus (strain ATCC 76670 / Carlsberg bottom yeast no.2 / CBS 1503 / CLIB 180 / NBRC 10610 / NRRL Y-1525) (Saaz-type lager yeast) (Saccharomyces monacensis) ORM1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces pastorianus (strain ATCC 76670 / Carlsberg bottom yeast no.2 / CBS 1503 / CLIB 180 / NBRC 10610 / NRRL Y-1525) (Saaz-type lager yeast) (Saccharomyces monacensis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces pastorianus (strain ATCC 76670 / Carlsberg bottom yeast no.2 / CBS 1503 / CLIB 180 / NBRC 10610 / NRRL Y-1525) (Saaz-type lager yeast) (Saccharomyces monacensis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ORM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.