



OCA5 Antibody

Product Code	CSB-PA466569XA01SVH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B5VJQ0
Immunogen	Recombinant Saccharomyces cerevisiae (strain AWRI1631) (Baker's yeast) OCA5 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain AWRI1631) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain AWRI1631) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	OCA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.