



DAO1 Antibody

Product Code	CSB-PA006494XA01GHM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99042
Immunogen	Recombinant Trigonopsis variabilis (Yeast) DAO1 protein
Raised In	Rabbit
Species Reactivity	Trigonopsis variabilis (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Trigonopsis variabilis (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	DAO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.