



HIS3 Antibody

Product Code	CSB-PA146965XA01FZS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02986
Immunogen	Recombinant Lachancea kluyveri (strain ATCC 58438 / CBS 3082 / CCRC 21498 / NBRC 1685 / JCM 7257 / NCYC 543 / NRRL Y-12651) (Yeast) (Saccharomyces kluyveri) HIS3 protein
Raised In	Rabbit
Species Reactivity	Lachancea kluyveri (strain ATCC 58438 / CBS 3082 / CCRC 21498 / NBRC 1685 / JCM 7257 / NCYC 543 / NRRL Y-12651) (Yeast) (Saccharomyces kluyveri)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lachancea kluyveri (strain ATCC 58438 / CBS 3082 / CCRC 21498 / NBRC 1685 / JCM 7257 / NCYC 543 / NRRL Y-12651) (Yeast) (Saccharomyces kluyveri)
Lead Time	Made-to-order (14-16 weeks)
Target Names	HIS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.