



HIS5 Antibody

Product Code	CSB-PA346706XA01CZJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56099
Immunogen	Recombinant Candida maltosa (Yeast) HIS5 protein
Raised In	Rabbit
Species Reactivity	Candida maltosa (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida maltosa (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	HIS5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.