




HSH49 Antibody

Product Code	CSB-PA859497XA01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99181
Immunogen	Recombinant Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) HSH49 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Target Names	HSH49

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: HSH49 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 28.6 kDa

Observed band size: 29 kDa