



IDH1 Antibody

Product Code	CSB-PA517526XA01AZC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13302
Immunogen	Recombinant Ajellomyces capsulatus (Darling's disease fungus) (Histoplasma capsulatum) IDH1 protein
Raised In	Rabbit
Species Reactivity	Ajellomyces capsulatus (Darling's disease fungus) (Histoplasma capsulatum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ajellomyces capsulatus (Darling's disease fungus) (Histoplasma capsulatum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IDH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.