



ADK2 Antibody

Product Code	CSB-PA269080XA01SAC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E9P974
Immunogen	Recombinant Saccharomyces cerevisiae (Baker's yeast) ADK2 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ADK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.