



uppS Antibody

Product Code	CSB-PA322605XA01SOE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20182
Immunogen	Recombinant Streptomyces fradiae (Streptomyces roseoflavus) uppS protein
Raised In	Rabbit
Species Reactivity	Streptomyces fradiae (Streptomyces roseoflavus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptomyces fradiae (Streptomyces roseoflavus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	uppS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.