



LAC1 Antibody

Product Code	CSB-PA741103XA01GHN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70KY3
Immunogen	Recombinant Melanocarpus albomyces LAC1 protein
Raised In	Rabbit
Species Reactivity	Melanocarpus albomyces
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Melanocarpus albomyces
Lead Time	Made-to-order (14-16 weeks)
Gene Names	LAC1