



KRR1 Antibody

Product Code	CSB-PA516785XA01SVR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E7QBZ1
Immunogen	Recombinant Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast) KRR1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Saccharomyces cerevisiae (strain Zymaflore VL3) (Baker's yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	KRR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.