



206 Antibody

Product Code	CSB-PA926706XA01GHO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q853W0
Immunogen	Recombinant Mycobacterium phage Omega (Mycobacteriophage Omega) 206 protein
Raised In	Rabbit
Species Reactivity	Mycobacterium phage Omega (Mycobacteriophage Omega)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mycobacterium phage Omega (Mycobacteriophage Omega)
Lead Time	Made-to-order (14-16 weeks)
Target Names	206
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.