



LIP5 Antibody

Product Code	CSB-PA868468XA01CZD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P8W0
Immunogen	Recombinant Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast) LIP5 protein
Raised In	Rabbit
Species Reactivity	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Lead Time	Made-to-order (14-16 weeks)
Target Names	LIP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.