



Malate synthase, glyoxysomal Antibody

Product Code	CSB-PA325532XA01RMM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17815
Immunogen	Recombinant Ricinus communis (Castor bean) Malate synthase, glyoxysomal protein
Raised In	Rabbit
Species Reactivity	Ricinus communis (Castor bean)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ricinus communis (Castor bean)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.