



# MPS1 Antibody

|                            |                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA848339XA01EKH                                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                              |
| <b>Uniprot No.</b>         | Q8SSH4                                                                                     |
| <b>Immunogen</b>           | Recombinant Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite) MPS1 protein |
| <b>Raised In</b>           | Rabbit                                                                                     |
| <b>Species Reactivity</b>  | Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)                          |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                               |
| <b>Form</b>                | Liquid                                                                                     |
| <b>Conjugate</b>           | Non-conjugated                                                                             |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4           |
| <b>Purification Method</b> | Antigen Affinity Purified                                                                  |
| <b>Isotype</b>             | IgG                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                 |
| <b>Product Type</b>        | Polyclonal Antibody                                                                        |
| <b>Immunogen Species</b>   | Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)                          |
| <b>Lead Time</b>           | Made-to-order (14-16 weeks)                                                                |
| <b>Gene Names</b>          | MPS1                                                                                       |